

## **Alltech® RAPIREAD™**

### **Customer User Manual**



#### **What is Alltech® RAPIREAD™?**

Alltech® RAPIREAD™ is the latest addition to our integrated system of tools and technologies that provide the most comprehensive protection from the threat of mycotoxins. Suitable for testing individual grains and forages, the system highlights real-time mycotoxin contamination on-farm or in the feed mill, providing you with 24/7 access to detailed results and analysis.

#### **The Alltech RAPIREAD system consists of:**

- The Raptor® Integrated Analysis Platform.
- Alltech RAPIREAD Mobile App.
- Alltech RAPIREAD Data Sync Tool.
- Mycotoxin Management Portal.

Accurate results can currently be determined in the field or the feed mill environment for aflatoxin, ochratoxin, zearalenone, DON, T2/HT-2 toxins, fumonisin and ergot toxins.

**Important: Before using Alltech RAPIREAD, please create an Alltech Mycotoxin Management Portal account by following the steps in this separate [guide](#). Once your account is set up, you can use the instructions in this manual to proceed to submit a test request and carry out mycotoxin testing with Alltech RAPIREAD.**

## **Table of Contents**

**Page 1.** Downloading and installing the Data Sync Tool

**Page 4.** Downloading and installing the Alltech RAPIREAD mobile app

**Page 7.** Carrying out a mycotoxin test (Submitting the test request)

**Page 9.** Running the mycotoxin test on the Neogen Raptor test device

**Page 11.** Uploading test results using the Alltech RAPIREAD Data Sync Tool

**Page 16.** Viewing mycotoxin quick reports on the Alltech RAPIREAD mobile App

**Page 19.** Other features of the Mycotoxin Management Portal


- Navigating the portal
- My Data Dashboard and reports
- Geographic overview

**Page 27.** Changing settings and managing users


**Page 30.** Alltech RAPIREAD support

## Step 1: Download and install the Alltech RAPIREAD Data Sync Tool on your PC or laptop


- Visit the Mycotoxin Management Portal at [www.knowmycotoxins.com](http://www.knowmycotoxins.com) and click CLIENT LOGIN.
- Enter your user credentials and click Login (this is username and password you will have used for setting up your portal account).
- When you arrive at the home page, click **Alltech RAPIREAD**, found on the orange header, and then select **Data Sync Tool**. Wait for the setup file named 'publish.zip' to download.



- Once downloaded, click the download icon in the top left corner. If a Google Drive warning appears during the downloading process, click **download anyway**.



- E. Once the file has downloaded, open the publish folder and double click the setup file. If a Windows Defender warning appears, click **More info** and click **Run anyway**.




- F. Click **Install** and wait until the program fully installs.
- G. The Alltech RAPIREAD Data Sync Tool should now be installed, and you will be taken to the sign-up screen.
- H. You will now be able to access the Alltech RAPIREAD Data Sync Tool from your windows desktop under **Alltech RAPIREAD Data Sync Tool**.

**Signing into the Alltech RAPIREAD Data Sync Tool:**


1. Open the Alltech RAPIREAD Data Sync Tool on your PC.
2. Enter your username and password (These are the same credentials as the Mycotoxin Management Portal).
3. Click Submit.
4. Go to your emails and wait for the verification code email. Copy and paste the verification code into the appropriate fields.
5. Click Sign In.
6. The Alltech RAPIREAD Data Sync Tool will now open and you can begin to sync files.

**Step 2: Downloading the Alltech RAPIREAD Mobile App on your iOS or Android mobile device**


- A. In the iOS App Store or Google Play store, search “Alltech RAPIREAD”.
  - a. For iPhone users, click GET in the iOS App Store and enter your Apple ID and username to install.
  - b. For Android users, click Install in the Google Play Store to install.




- B. Once installed, open the app and you will be taken to the sign in page. Sign in with your portal email and password.



C. A verification code will be sent to your email.



D. Enter the verification code and click **Verify**.



Congratulations! You have now downloaded and installed each of the necessary tools to run mycotoxin tests and view the subsequent results on Alltech RAPIREAD.

***Note: You will be required to sign in every 30 days. The previous login will expire after 30 days, and you will require a new verification code when signing in again.***


## Carrying out a mycotoxin test with Alltech RAPIREAD

Note: The same Mycotoxin Management Portal login credentials must be used for the Alltech RAPIREAD Data Sync Tool and the Alltech RAPIREAD Mobile App to successfully carry out a test and transfer the information.

### Step 1: Create a test request using the Alltech RAPIREAD Mobile App

- A. Open the Alltech RAPIREAD Mobile App and login if required.
- B. Click on the Test Request tab at the bottom of the screen.
- C. Complete all the required fields.

*Tip: All the fields on the Test Request tab can be prefilled to save time. For example: If a test is always going to be conducted in a certain State/Region, that box can be prefilled. You can also default to accept the Alltech RAPIREAD validation warnings. See the Company Settings on page 28 for defaulting Alltech RAPIREAD Test Request Fields on the Alltech RAPIREAD Mobile App.*

A screenshot of the Alltech RAPIREAD Mobile App interface. It shows a dark grey header bar on the left. The main content area is white and contains three checked checkboxes with the following text: "By default, I accept the Forages validation warning.", "By default, I accept the Ingredients validation warning.", and "By default, I accept the Mixed Feeds validation warning.". Below these checkboxes is a light blue horizontal bar. At the bottom of the form are two buttons: a green "Save" button and a grey "Cancel" button.

- D. Click Submit.
- E. Check that all details are correct on the submission page and click **Confirm** to finish.

Confirm Submission

ORIGIN OF SAMPLE (COUNTRY)  
United States

STATE/REGION  
Indiana

CITY  
Seymour

CUSTOMER SAMPLE ID (OPTIONAL)  
Corn Bin 1

CUSTOMER  
Office - Missouri

CONTACT  
Josh Teulker

SPECIES CATEGORY  
Poultry

SPECIES  
Layers

FEED CATEGORY  
Ingredients


FEED INGREDIENT  
Corn

REPORT LANGUAGE  
English

Confirm

- F. A QR Code will be generated. This code will be used on the Raptor in future steps.

Test Request Submission Form



SUBMISSION ID:  
29-28363437-336

New Test Request

Home Test Request Submissions

*Tip: If you lose the QR code you created, go to the **Submissions** tab on the Alltech RAPIREAD Mobile App and click on the Test Request. This will bring the QR code back up on the app.*


*Note: Ensure you have good internet or mobile connection when creating a test request. If connection is interrupted, the test request may not work properly.*

## Step 2: Carrying out the mycotoxin test on the Neogen Raptor device

- A. Follow the guidelines that Neogen provides to prepare a sample for testing. These guidelines are found in the booklets provided in each Reveal Quantitative Test Box.



- B. Once the sample is prepared and ready to test, insert the appropriate test strip in a cartridge and insert the cartridge into 1 of the 3 ports on the top of the Raptor.
- C. The Raptor will read the test strip and ask to either Enter Sample ID or Scan Sample ID. Choose **Scan Sample ID**.



- D. Scan the QR Code that was created on the Alltech RAPIREAD Mobile App for the sample being tested.



***Tip: Ensure the backlight on your phone is turned to the brightest setting.***

- E. Continue to follow Neogen's instructions to complete the mycotoxin test.

***Note: Each QR Code represents 1 sample. Continue to scan the same QR code if testing for multiple mycotoxins. You will require 3 or more toxin tests to get all features, such as REQ, on the quick report.***


### **Step 3: Uploading test results using the Alltech RAPIREAD Data Sync Tool**

- A. Once all mycotoxin tests have been completed for the sample, use the USB cable to connect the Neogen Raptor device to your laptop or PC.


*Note: This step will only work via USB. Network connection will result in an error.*




- B. Open and sign in to the Alltech RAPIREAD Data Sync Tool if necessary.  
C. On the Raptor, select the computer icon and select the **Connect with USB** option.



- D. It may take a few minutes for the Raptor to install. Wait until your PC displays the **Device is ready** message before continuing.
- E. Once installed, click the **Connect** button in the right corner of the screen. All files on the Raptor should now display under the **Files** tab.




- F. Select the files you want to sync by clicking on the check box on the left side of the file name.



**Note: An error will occur if any test is selected that has already been uploaded to Alltech RAPIREAD or does not have a QR Code associated with the sample.**

G. To check the details of a test sample, click the **Details** tab.


H. Click **Upload** and wait for the files to finish syncing.



I. Your report will now be available to view.

## Step 4: Viewing mycotoxin quick reports on the Mycotoxin Management Portal

- Visit [www.knowmycotoxins.com](http://www.knowmycotoxins.com).
- Click on **Client Log In** and enter your Username and Password.
- Click **Alltech RAPIREAD** on the task bar and click **Individual Results**.




- Click **View Reports** under the Operations column.


Date	Customer	Contact	Sample Submission I	Customer Sample ID	Internal ID	Status	Operations
2020-05-05	Alltech	Office Contact	29-87223739-336	Corn 1	Rapid Test	Completion	View Test Request   View Reports   Generate Report
2020-04-21	Alltech	Office Contact	29-83323469-336	Multiple test	Rapid Test	Completion	View Test Request   View Reports   Generate Report
2020-04-21	Alltech	Office Contact	29-80265370-336	Corn test 2	Rapid Test	Completion	View Test Request   View Reports   Generate Report
2020-04-17	Alltech	Office	29-28547852-336	Corn	Rapid Test	Completion	View Test Request   View Reports   Generate Report

- Scroll to the bottom of the page. Under Final Reports, click **Assessment Report** under the Reports column.

# Altech<sup>®</sup> MYCOTOXIN MANAGEMENT




*Note: You can also generate new reports for different species by clicking **Generate Report** under the operations column on the home screen. Select a new species in the dropdown box and click **Generate Report** in the bottom right corner. Your new report will appear under Final Reports.*




**Step 5: Viewing a Mycotoxin Quick Report on the Alltech RAPIREAD Mobile App**

- A. Open the Alltech RAPIREAD Mobile App and sign in if required.
- B. Click on the **Submissions** tab at the bottom of the screen.



- C. From here, you can view a full list of all your mycotoxin test submissions.




Verizon 12:44 PM

**Submissions**

DATE	SUBMISSION ID / CUSTOMER SAMPLE ID	STATUS	REPORT
01/14/2021	29-41453469-336 Corn	✓	View
01/08/2021	29-29403412-336 Corn bin 4a retest	✓	View
12/23/2020	29-36304183-336 12/20/20 Corn Bin 4B	✓	View
12/23/2020	29-35425780-336 12/20/20 Corn Bin 2-B	✓	View
12/23/2020	29-34500877-336 12/20/20	✓	View
12/23/2020	29-33004653-336 12/20/20 Corn Bin 3-A	✓	View
12/23/2020	29-32055334-336 12/20/20 Corn Bin 4-A	✓	View

Home Test Request Submissions

D. Under the report column, click **View** to view the results of a mycotoxin test.



Verizon 12:44 PM

**Submissions**

DATE	SUBMISSION ID / CUSTOMER SAMPLE ID	STATUS	REPORT
01/14/2021	29-41453469-336 Corn	✓	View
01/08/2021	29-29403412-336 Corn bin 4a retest	✓	View
12/23/2020	29-36304183-336 12/20/20 Corn Bin 4B	✓	View
12/23/2020	29-35425780-336 12/20/20 Corn Bin 2-B	✓	View
12/23/2020	29-34500877-336 12/20/20	✓	View
12/23/2020	29-33004653-336 12/20/20 Corn Bin 3-A	✓	View
12/23/2020	29-32055334-336 12/20/20 Corn Bin 4-A	✓	View

Home Test Request Submissions

*Note: This can only be done when the submission from the Alltech RAPIREAD Mobile App and the results from the Raptor have been synced. If the status column reads 'pending,' you will need to sync the Alltech RAPIREAD Mobile App and results.*

E. Click **Back** to navigate back to the view of all your test submissions.

Mycotoxin Quick Report

Customer:

Species: Dairy Cows

Sample Submission: 29-29403412-336

Matrix: Corn

Customer Sample ID: Corn bin 4a retest

Mycotoxins	Level, ppb	Lower	Moderate	Higher
Deoxynivalenol	16600	500	1000	2000

Mycotoxin Assessment for Performance Impairment

Deoxynivalenol

Estimated Risk\*

Higher

Potential Impact on Dairy Cow Performance

Milk production loss, lb/cow/day (1 lb/cow milk = 2,272 kcal)

0 -0.4 -0.5 -0.6 -0.9 -1.4

Change in SCC, %

0 10 20 80 160 725

Comments:

The Alltech Quick REQ (Risk Equivalent Quantity) of this feedstuff represents a Higher risk level to animal performance. The overall dietary risk will depend on the inclusion rate of this feedstuff and mycotoxin concentration of other raw materials. Mycotoxins may alter feed intake, gut health, immunity and productivity.

Generated 11/26/2021

alltech

MYCOTOXIN MANAGEMENT

Back

## **Other features of the Mycotoxin Management Portal**

The new Mycotoxin Management Portal offers many new functions to users, including a My Data Dashboard where users can filter samples by different species and generate results for different Data Summaries, Risk Assessments, Risk Timelines and Risk Maps. More features include the ability to view a Geographic Overview of the samples tested and export data to an excel sheet.

## **Navigating the new Mycotoxin Management Portal**

The Home screen now gives an overview of all Mycotoxin tests conducted with Alltech RAPIREAD and Alltech 37+® over the past three months, as well as allowing you to switch to any of the six available languages on the portal. From the Home screen, you can navigate to the new features of the Mycotoxin Management Portal, including My Data Dashboard, Geographical Overview and Exporting Data. For instructions to view Individual mycotoxin test results, see step 4 of the manual.

## **Viewing My Data Dashboard**

- A. From the Home screen, click the Alltech RAPIREAD button on the orange ribbon.
- B. Then click My Data Dashboard
- C. On the My Data Dashboard screen, select your language and a date range you would like to view samples in by entering the dates in the gray boxes or by moving the orange circles on the slider bar.
- D. You can also filter between Customer ID, Sample ID, Ingredient/Complete feed, Country and State/Province/County if you have a specific need for viewing specific samples.
- E. Next, select the species for the risk assessment at the bottom of the page.

Alitech MYCOTOXIN MANAGEMENT

United States - English

English Español Français Português Русский 中文

**Select samples**

1/1/2021 4/22/2021

Sample ID Country

All All

Customer ID Ingredient/Complete feed State/Province/County

All All All

**Select species for risk assessment**

Pigs - Growers/Finishers

**Go to results**

Data summary


Risk assessment

Risk timelines


Risk map

F. Once all selections for the samples are complete, you can then click any of the boxes on the right side of your screen to go to your results.


G. There are four options to view results: Data Summary, Risk Assessment, Risk Timelines and a Risk Map. To view different pages, simply click on the different options on the left side of the page. A description of each view is shown below.



**Data Summary:** Shows the number of samples analyzed and mycotoxins tested, as well as which toxins were tested. There is also a bar chart showing the average mycotoxin concentration for the selected samples or the maximums for each mycotoxin tested. To switch from Averages to Maximums, simply click on the correlating buttons to the left of the bar chart.




**Risk Assessment:** The Risk Assessment allows you to select the ingredient inclusion rate for the estimated risk to the ration for the samples selected. Simply complete all the drop-down boxes to the left of the chart to view the REQ of samples tested.




\*Values shown may be below the limit of quantification (LOQ) for the test kit. Blank values indicate mycotoxin not tested. Risk assessment based on Alltech Practical Limits.

**Mycotoxin Risk Timeline:** Allows the user to view the risk of each mycotoxin tested over the selected period of time. Simply select the mycotoxin you would like to see the risk for and the species you are feeding to view the risk level for each month. The chart will default to show Alltech's upper and lower practical limits for the selected species. To view the chart without the limits, simply click No under Show Alltech practical limits?




**Mycotoxin Risk Maps:** This allows you to view your mycotoxin results via an interactive map showing the risk level of samples by location, based on the location entered on the submission form. Simply select the species you are feeding from the dropdown to the left of the map and click which toxin you would like to view on the top of the map.



**Viewing Geographic Overview:** This view allows users to see multiple function on one screen, including an interactive map, REQ bar chart for the different mycotoxins tested, the number of samples analyzed and mycotoxin tests performed for the calendar year selected, giving users a summary of all samples submitted by every user during the time period selected. These samples are displayed anonymously to ensure sensitive data is not available to all users.

- A. On the Home screen, click the Alltech RAPIREAD tab then click Geographical Overview.
- B. On the Geographical Summary screen, select your language, then change any of the settings on the left side of the screen to view a different calendar year, species fed, ingredient and mycotoxin tested.



## Exporting User Data to an Excel file (Users can only export their own test data)

- A. On the Home screen, click the Alltech RAPIREAD Tab then click Export Data.
- B. Enter in the date range you would like to export and select the Feed Matrix.
- C. Click Download

Home Alltech 37+ - Alltech RAPIREAD -

Date: From 5/28/2020 Date: To 5/28/2021

Past 1 Month Past 3 Month Past 1 Year

Feed Matrix  
Barley


Download

- D. The Excel file will begin to download. Once downloading is complete, you can open the file to view your results.

TestSampleData (1) - Protected View																											
RequestID	SampleID	SubmitMatrix	Country	StateProv	CityRegion	Species	DateSampleReceived	Aflatoxins	Aflatoxins	Aflatoxins	Ochratoxin	Ochratoxin	Ochratoxin	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol	Deoxynivalenol
1	4001	2018 Barle Farm PL	Barley		Stavelly		12/9/2020	3.9						88													
2	4002	2019 Barle Farm PL	Barley		Stavelly		12/9/2020	2.4						113.9													
3	4003	Barley Farm TN	Barley		Ponoka		12/10/2020	2.4						19.2													
4	4004	Barley Farm RE	Barley		Olds		1/14/2021							52.3													
5	4005	Barley G01 Granum V	Barley		Granum		11/22/2020	0.99						32.5													
6	4006	Barley G01 Granum V	Barley		Granum		11/22/2020	1.1						42.6													
7	4007	Barley G01 Granum V	Barley		Granum		12/13/2020							186.5													
8	4008	Barley G01 Granum V	Barley		Granum		12/16/2020							40.1													
9	4009	Barley G01 Granum V	Barley		Granum		12/16/2020							161.8													
10	4010	Barley G01 Granum V	Barley		Granum		12/13/2020							151.8													
11	4011	Barley G01 Granum V	Barley		Granum		12/20/2020							131.6													
12	4012	Barley G01 Granum V	Barley		Granum		12/16/2020							133.1													
13	4013	Barley G01 Granum V	Barley		Granum		1/2/2021							89.3													
14	4014	Barley G01 Granum V	Barley		Granum		1/2/2021							101.7													
15	4015	Barley G01 Granum V	Barley		Granum		1/2/2021							106.5													
16	4016	Barley G01 Granum V	Barley		Granum		12/20/2020							196.7													
17	4017	Barley G01 Granum V	Barley		Granum		12/20/2020							132.5													
18	4018	Barley G01 Granum V	Barley		Granum		1/2/2021							182.7													
19	4019	Barley G01 Granum V	Barley		Granum		1/2/2021							214.1													
20	4020	Barley G01 Granum V	Barley		Granum		1/2/2021							76													
21	4021	Barley G01 Granum V	Barley		Granum		12/13/2020							137													
22																											
23																											
24																											
25																											
26																											
27																											

## Settings

The Settings tab allows you to Manage Users, Change Company Settings, Edit Account Settings and Change Password. Each function is outlined below.



**Manage Users:** Allows you to enable and manage the contacts that are associated with the account. To enable or disable an account, simply click on the red X under the enable column to disable or click on the green check to enable. The Representative of the account can enable any account associated with the portal. Contact [jteulker@alltech.com](mailto:jteulker@alltech.com) to have a Representative assigned to the account if the current Representative cannot be contacted.

**Alltech MYCOTOXIN MANAGEMENT**

United States - English

Home Alltech 37++ Alltech RAPIREAD - Support Settings Logout

Contacts

Search

Name	Username	Email	Phone	Enable	Representative
Maddy Lazurko	MTJF	mtjg-1154@alltech.com		✓	<input checked="" type="checkbox"/>
Maddy Lazurko	MTFJ	mtjg-1153@alltech.com	3060617487	✓	<input type="checkbox"/>

**Company Settings:** Allows you to make changes to your Name, Address and Phone Number. Under the company setting, you can also default any of your Test Request settings for the Alltech RAPIREAD Mobile App. For example, if a user is always feeding the same species, default the Species Category and Species to that species.

### Defaulting the Alltech RAPIREAD Mobile App Test Request

- A. Click on the Settings tab in the top right then click Company Settings.
- B. On the Company Settings page, click Alltech RAPIREAD in the bottom left corner.

**Alltech MYCOTOXIN MANAGEMENT**

Home Alltech 37+ Alltech RAPIREAD -

Customer Name	Country
MTJF	Canada
City	Address1
Strathmore	126 Canal Gardens
Phone	Fax
18669484578	
Alltech RAPIREAD	
Save	

- C. In the Alltech RAPIREAD popup, you can then default any of the Test Request form boxes. Ensure you also check the default boxes on the Validation Warnings to prevent the warning from appearing for each test request performed. Once you have filled in the boxes, read and accept the terms and conditions and click Save.

The screenshot shows a web form titled "Alltech RAPIREAD" with a close button (X) in the top right corner. The form is set against a background with tabs for "Canada" and "Alberta". The form contains the following fields:

- Customer Sample ID (optional): Text input field.
- Language: Dropdown menu with "English" selected.
- Country: Text input field with "Canada" entered.
- State: Dropdown menu with "Alberta" selected.
- City: Text input field with "Strathmore" entered.
- Species Category: Dropdown menu.
- Species: Dropdown menu.
- Species Other Txt: Text input field.
- Feed Category: Dropdown menu.
- Feed: Dropdown menu.
- Feed Other Txt: Text input field.

At the bottom left, three checkboxes are grouped together and circled in red:

- ☒ By default, I accept the Forages validation warning.
- ☒ By default, I accept the Ingredients validation warning.
- ☒ By default, I accept the Mixed Feeds validation warning.


At the bottom of the form are two buttons: "Save" (green) and "Cancel" (grey).

**Account Settings:** Allows users to change the name and phone number associated with their account.

**Change Password:** Allows users to change the password associated with their account.

## Alltech RAPIREAD support

For answers to any questions regarding the Mycotoxin Management Portal, please use the Support tab on the top right side of the screen. There you can find a list of Frequently Asked Questions or a link to Contact Us directly.



## Instructional video and further questions

You can view a short video at this [link](#) that outlines the main steps required to submit a mycotoxin test, run the test on the Raptor device and view the results on your portal and Alltech RAPIREAD Mobile App.

If you have further questions on getting set up or operating Alltech RAPIREAD, please contact [mycotoxinsupport@alltech.com](mailto:mycotoxinsupport@alltech.com)