







Reveal Q+ MAX Products

<p>Industry leading tests for mycotoxins</p> <ul style="list-style-type: none"> • Water based extraction • A single extraction for all 6 mycotoxins • Results in minutes 	Aflatoxin	Ochratoxin	Zearalenone	T2/HT-2	Fumonisin	DON
						
Part Number	8088	8688	8188	8288	8888	8388
Limit of Detection	2 ppb	1.10 ppb	corn: 36 ppb wheat:21 ppb	40 ppb	0.25 ppb	0.3 ppm
Range of Quantitation (with dilutions)	2-300 ppb	2-100 ppb	corn: 25-1500 ppb wheat: 50-1500 ppb	50-3000 ppb	0.25-50 ppm	0.3-30 ppm
Time to result	6 min	5 min	5 min	5 min	5 min	3 min
GIPSA	Approved	being submitted	being submitted	NA	being submitted	being submitted
Key Features	best cross reactivity on the market B ₁ 100% B ₂ 91% G ₁ 87% G ₂ 64%	first LF validated on green coffee fastest time to result of LF kits on market	first LF with aqueous extraction fastest testing time	validated on several grain commodities - competitors have only corn	common extraction will allow users to run several mycotoxins off one cup - saves time and money on testing supplies	common extraction will allow users to run several mycotoxins off one cup - saves time and money on testing supplies

*common extraction - each toxin will require a different diluent so you do have to purchase each toxin kit separately

*Kit includes a box of 25 extraction packets - you may want to consider purchasing an extra box of MAX1 Extraction packets in case you have to repeat an assay. These are available as Neogen part # 8089.

*validation reports are not currently available but coming soon

*support literature for each of the available kits is located on Tech Data - will be adding additional information as the kits become available.

*Kits are for use with AccuScan GOLD and PRO readers only at this time